



ZYMO RESEARCH

DNA
Purification
MADE SIMPLE
Made Simple™

Quick-DNA™ Fungal/Bacterial Miniprep Kit

DNA from tough-to-lyse fungi and bacteria samples.

Highlights

- Simple, efficient isolation of DNA (up to 25 µg/prep) from all types of tough-to-lyse fungi (e.g., yeast) and bacteria in as little as 15 minutes.
- State-of-the-art, ultra-high density **BashingBeads™** are fracture resistant and chemically inert.
- Omits the use of organic denaturants as well as proteinases.

Catalog Numbers:
D6005



Scan with your smart-phone camera to view the online protocol/video.



tech@zymoresearch.com



www.zymoresearch.com



Toll Free: (888) 882-9682

Table of Contents

Product Contents	01
Specifications	02
Product Description	03
Protocol	04
Ordering Information	06
Guarantee	08

Product Contents

Quick-DNA™ Fungal/Bacterial Miniprep Kit	D6005 (50 Preps.)	Storage Temperature
ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm)	50	Room Temp.
BashingBead™ Buffer	40 ml	Room Temp.
Genomic Lysis Buffer ¹	100 ml	Room Temp.
DNA Pre-Wash Buffer ²	15 ml	Room Temp.
g-DNA Wash Buffer	50 ml	Room Temp.
DNA Elution Buffer	10 ml	Room Temp.
Zymo-Spin™ III-F Filters	50	Room Temp.
Zymo-Spin™ IICR Columns	50	Room Temp.
Collection Tubes	150	Room Temp.
Instruction Manual	1	-

¹ For optimal performance, add beta-mercaptoethanol to 0.5%(v/v) *i.e.*, 500 µl per 100 ml.



² A precipitate may have formed in the **DNA Pre-Wash Buffer** during shipping. To completely resuspend the buffer, incubate the bottle at 30-37°C for 30 minutes and mix by inversion. DO NOT MICROWAVE.

Specifications

- **Format** – Bead Beating, Spin Column Purification
- **Sample Sources** – 50 - 100 mg (wet weight) fungi or bacteria; this equates to approximately 10^9 bacterial cells and 10^8 yeast cells. Spores, pollen, nematodes, as well as other microorganisms can also be sampled.
- **DNA Purity** – High quality DNA is eluted with **DNA Elution Buffer** making it perfect for PCR ($A_{260}/A_{280} >1.8$).
- **DNA Size Limits** – Capable of recovering genomic DNA up to and above 40 kb. In most instances, mitochondrial DNA and viral DNA (if present) will also be recovered.
- **DNA Recovery** – Typically, up to 25 μg total DNA is eluted into 100 μl (35 μl minimum) **DNA Elution Buffer** per sample.
- **Equipment** – Microcentrifuge, Vortex, Cell Disrupter/Pulverizer (recommended)

Product Description

The **Quick-DNA™ Fungal/Bacterial Miniprep Kit** is designed for the simple, rapid isolation of DNA from tough-to-lyse fungi, including *A. fumigatus*, *C. albicans*, *N. crassa*, *S. cerevisiae*, *S. pombe*, as well as from mycelium and Gram positive and Gram negative bacteria. The procedure is easy and can be completed in as little as 15 minutes: fungal and/or bacterial samples are added directly to a **ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm)** and rapidly and efficiently lysed by bead beating without using organic denaturants or proteinases. The DNA is then isolated and purified using our Zymo-Spin™ Technology and is ideal for downstream molecular-based applications including PCR, array, etc. A schematic of the **Quick-DNA™ Fungal/Bacterial Miniprep Kit** procedure is shown below.



DNA isolated from *Saccharomyces cerevisiae* (spores) and *E. coli* using the **Quick-DNA™ Fungal/Bacteria Kit** was high-quality and structurally intact. Equivalent amounts of yeast and bacteria were processed using the **Quick-DNA™ Fungal/Bacterial Kit** or the Supplier A kit. Equal volumes of eluted DNA were analyzed on a 0.8% (w/v) agarose gel stained with EtBr.

DNA/RNA Shield™ (R1100-50, R1100-250) can be used to stabilize nucleic acids and inactivate infectious agents in a variety of samples, without the need for reagent removal.

For rapid, robust, and simple purification of high quality, inhibitor-free DNA from any sample including feces, soil, water, biofilms, swabs, saliva, body fluids, etc. use the **ZymoBIOMICS™ DNA Miniprep Kit (D4300)**.

Protocol

For optimal performance, add beta-mercaptoethanol (user supplied) to the **Genomic Lysis Buffer** to a final dilution of 0.5 %(v/v) *i.e.*, 500 µl per 100 ml.

1. Add 50 – 100 mg (wet weight) fungal or bacterial cells¹ that have been resuspended in up to 200 µl of water or isotonic buffer (*e.g.*, PBS) to a **ZR BashingBead™ Lysis Tube (0.1 mm & 0.5 mm)**. Add 750 µl **BashingBead™ Buffer** to the tube².
2. Secure in a bead beater fitted with a 2 ml tube holder assembly and process at maximum speed for ≥ 5 minutes.

Note: Required processing time will vary depending on the device and application and therefore should be evaluated on a case by case basis.

*For example, processing times may be as little as 3 minutes when using high-speed cell disrupters (*e.g.*, the portable TerraLyzer™ Sample Processor, FastPrep® -24, or similar) or as long as 20 minutes when using lower speeds (*e.g.*, Disruptor Genie™, or standard benchtop vortexes). See manufacturer's literature for operating information.*

3. Centrifuge the **ZR BashingBead™ Lysis Tube (0.1 & 0.5 mm)** in a microcentrifuge at 10,000 x g for 1 minute.
4. Transfer up to 400 µl supernatant to a **Zymo-Spin™ III-F Filter** in a **Collection Tube** and centrifuge at 8,000 x g for 1 minute.
5. Add 1,200 µl of **Genomic Lysis Buffer** to the filtrate in the **Collection Tube** from Step 4.
6. Transfer 800 µl of the mixture from Step 5 to a **Zymo-Spin™ IICR Column³** in a **Collection Tube** and centrifuge at 10,000 x g for 1 minute.
7. Discard the flow through from the **Collection Tube** and repeat Step 6.
8. Add 200 µl **DNA Pre-Wash Buffer** to the **Zymo-Spin™ IICR Column** in a new **Collection Tube** and centrifuge at 10,000 x g for 1 minute.
9. Add 500 µl **g-DNA Wash Buffer** to the **Zymo-Spin™ IICR Column** and centrifuge at 10,000 x g for 1 minute.

¹ This equates to approximately 10⁹ bacterial cells and 10⁶ yeast cells.

² Cap tube tightly to prevent leakage.

³ The **Zymo-Spin™ IICR Column** has a maximum capacity of 800 µl.

10. Transfer the **Zymo-Spin™ IICR Column** to a clean 1.5 ml microcentrifuge tube and add 100 μl (35 μl minimum) **DNA Elution Buffer** directly to the column matrix. Centrifuge at 10,000 x *g* for 30 seconds to elute the DNA.

Ultra-pure DNA is now ready for use in your experiments.

Ordering Information

Product Description	Catalog No.	Size
Quick-DNA™ Fungal/Bacterial Microprep Kit	D6007	50 Preps.
Quick-DNA™ Fungal/Bacterial Miniprep Kit	D6005	50 Preps.
Quick-DNA™ Fungal/Bacterial Midiprep Kit	D6105	25 Preps.
Quick-DNA™ Fungal/Bacterial 96 Kit	D6006	2 x 96 Preps.

Individual Kit Components	Catalog No.	Amount
Genomic Lysis Buffer	D3004-1-100	100 ml
BashingBead™ Buffer	D6001-3-40	40 ml
DNA Pre-Wash Buffer	D3004-5-15	15 ml
g-DNA Wash Buffer	D3004-2-50	50 ml
DNA Elution Buffer	D3004-4-10	10 ml
ZR BashingBead™ Lysis Tubes (0.1 & 0.5 mm)	S6012-50	50 Tubes
Zymo-Spin™ III-F Filters	C1057-50	50 Pack
Zymo-Spin™ IICR Columns	C1078-50	50 Pack
Collection Tubes	C1001-50 C1001-500 C1001-1000	50 Pack 500 Pack 1,000 Pack

Lysis Instruments	Catalog No.	Amount
-------------------	-------------	--------



Disruptor Genie™, 120V w/ 2 ml tube holder assembly.	S6001-2-120	1 Unit
Disruptor Genie™, 230V w/ 2 ml tube holder assembly.	S6001-2-230	1 Unit
TurboMix Attachment, 2 ml Permanently mounts to most existing Vortex Genie™ mixers converting them to a Disruptor Genie™.	S6004	1 Unit

The **Disruptor Genie™** with 2 ml tube holder assembly from Scientific Industries, Inc. (Cat. No. S6001-2-120 from Zymo Research Corp.)



**100% satisfaction guarantee on all Zymo Research products,
or your money back.**

Zymo Research is committed to simplifying your research with quality products and services. If you are dissatisfied with this product for any reason, please call 1(888) 882-9682.

Integrity of kit components is guaranteed for up to one year from date of purchase.
Reagents are routinely tested on a lot-to-lot basis to ensure they provide the highest performance and reliability.

This product is for research use only and should only be used by trained professionals. It is not for use in diagnostic procedures. Some reagents included with this kit are irritants. Wear protective gloves and eye protection. Follow the safety guidelines and rules enacted by your research institution or facility.

[™] Trademarks of Zymo Research Corporation
Disruptor Genie[™] is a trademark of Scientific Industries, Inc. and FastPrep[®] is a registered trademark of Qbiogene, Inc.



ZYMO RESEARCH

The **BEAUTY** of **SCIENCE** is to Make Things **SIMPLE**®



tech@zymoresearch.com



www.zymoresearch.com



Toll Free: (888) 882-9682