



Vicbio Biotechnology Co., Ltd.
Block A, No.2, Jingyuan Street, Beijing, China, 100176
www.nebulabio.com

Certificate of Analysis

Product Name: N1-Methylpseudo-UTP Sodium Salt, 100mM
Catalog Number: NBM007
Lot Number: BA10000
Date of Manufacture: Oct.20.2021
Country of Manufacture: China

Test	Specification	Result
1. Appearance	Clear, colorless liquid	Pass
2. pH	7.0±0.1	7.1
3. Purity (HPLC area%)	≥99%	99.7%
4. Concentration	100mM±5%	103mM
5. DNase	0	Pass
6. RNase	0	Pass
7. Quality Control	Molecular Biology Grade This lot of NTP has been performance tests in transcription reaction.	Pass
8. Storage	-20°C. Avoid repeated freeze-thaw cycles. Upon first use, it is recommended to aliquot the reagent into single use portions.	N/A

Restest Date: Nov, 2022