


Product Datasheet

Product Number **NBM008**
Revision Number **NV001**



Product Name	5-Methoxy-UTP Sodium Salt, 100mM
Alias	5moUTP; 5-Methoxyuridine-5'-triphosphate, Sodium Salt
Molecule Formulation	C10H14N2Na3O16P3
Molecule Weight	580
CAS#	N/A
Purity	≥ 97 % (HPLC)
Form	Solution in water
Color	Colorless
Concentration	100mM±0.5%
pH (25°C)	7.0±0.1
Spectroscopic Properties	λ _{max} 279nm±2nm (pH 7.0)