


# Product Datasheet

Product Number **NBM001**  
Revision Number **NV001**



Product Name	ARCA Cap Analog Sodium Salt, 100mM
Alias	3'-O-Me-m7G(5')ppp(5')G; N7-Methyl-3'-O-Methyl-Guanosine-5'-Triphosphate-5'-Guanosine; ARCA Analog
Molecule Formulation	C <sub>22</sub> H <sub>32</sub> N <sub>10</sub> O <sub>18</sub> P <sub>3</sub> (free acid)
Molecule Weight	817.4 (free acid)
CAS#	N/A
Purity	≥ 95 % (HPLC)
Form	Solution in water
Color	Colorless
Concentration	100mM±4%
Storage	At or below -20°C