



BioWhittaker™ Specialty MediaMaximize Your Cell Growth and Productivity



Maximize Your Cell Growth and Productivity

If your work leads to downstream therapeutic, diagnostic, or biopharmaceutical applications, serum-free media are the best choice for initiating your research studies.

BioWhittaker™ Specialty Media deliver outstanding performance, a broad selection, and consistent results, whether you need proteinfree, non-animal origin (NAO) or chemically defined media. Our ISO 9000 and ISO 13485 systems are recognized industry wide.

BioWhittaker™ Serum-free Media deliver unique benefits:

- Consistent results
- Elimination of FBS lot qualification
- Elimination of FBS-borne mycoplasma and virus contamination
- Simplified downstream purification
- Reduced regulatory burdens
- Maximal yields
- Produced according to cGMP

General Purpose Media

UltraCULTURE™ Serum-free Medium

- Delivers excellent growth and expansion of a broad range of adherent and suspension cells
- Simplify your media needs, no FBS needed
- The formulation for UltraCULTURE™ Media has been submitted to the FDA as a Product Master File

PC-1™ Chemically Defined, Serum-free Medium

- Superior medium for primary adherent cells

HL-1™ Chemically Defined, Serum-free Medium

- Unique medium for hybridoma and lymphocytic cells

	General use	Non-animal origin	Protein-free	Chemically defined	Optimized for
UltraCULTURE™ Medium		-	-	-	Most cell types
PC-1™ CD Medium		_			Adherent primary cells and cell lines
HL-1™ CD Medium		_			Hybridomas and many adherent cell lines

Ordering Information

Cat. No US	Cat. No EU	Product	Size
BP12-725F	BEBP12-725F	UltraCULTURE™ Medium	500 mL
BP77232	BP77232	PC-1™ CD Medium	2 X 500 mL
BP77201	BP77201	HL-1™ CD Medium	2 X 500 mL

Mesenchymal Stem Cell Medium

MSCGM™, Serum-free MSC Growth Medium

- Optimized for multiple passage expansion of all types of hMSCs
- Cells can be directly transitioned from serum-containing medium with little to no adaptation time
- No need for attachment matrix to plate cells
- Supports multi-lineage differentiation

Ordering Information

Cat. No.	Product	Size
190632	TheraPEAK™ MSCGM™ Mesenchymal Stem Cell — Growth Medium BulletKit™	Kit
190620	TheraPEAK™ MSCBM™ Mesenchymal Stem Cell — Basal Medium	500 mL
192125	TheraPEAK™ MSCGM™ Mesenchymal Stem Cell — Growth Medium SingleQuots™ Kit	5 mL

Cryopreservation Medium

ProFreeze™ Chemically Defined Medium (2X)

- Non-animal, protein-free freezing medium
- Chemically defined

- Formulated for cryopreserving all cell types
- Maintains high cell viability upon recovery from frozen storage
- Requires 15% DMS0 at 2X

Ordering information

Cat. No US	Cat.No EU	Product	Size
BP12-769E	BP12-769E	ProFreeze™ Chemically Defined Medium (2X)	100 mL

Trypsinisation

TrypZean™ Trypsin/EDTA (10x)

- Non-animal origin
- Recombinant bovine trypsin expressed in corn

Ordering information

Cat. No US	Cat. No EU	Product	Size
BE02-034E	BE02-034E	TrypZean™ Trypsin/EDTA (10x) research use only (RUO)	100 mL

Hematopoietic Media

X-VIVO™ Xeno-free, Serum-free Hematopoietic Cell Media

- Superior media for most hematopoietic cell types including dendritic cells, lymphocytes, monocytes, macrophages, and granulocytes
- Some formulations include recombinant transferrin to meet the highest regulatory requirements
- All current TheraPEAK™ X-VIVO™ media products are manufactured under current GMPs and are listed with the FDA in a Product Master File.
- Highly cited in scientific literature

	General use	Non-animal origin	Protein-free	Optimized for
X-VIV0™ 10 Medium	-	-	-	PBL, LAK, monocytes, macrophages, stem cells
X-VIV0™ 15 Medium			_	PBL, NK, TIL, monocytes, macrophages, stem cells, dendritic cells
X-VIV0™ 20 Medium	_	_	_	PBL, LAK, TIL, monocytes, macrophages, stem cells, dendritic cells

Ordering information

Cat. No US	Cat. No EU	Product	Size		
04-3800	BE 04-380Q	X-VIVO™ 10 Medium With L-glutamine, Gentamicin and Phenol Red, With Human Transferrin	1 L		
BP04-743Q	BEBP04-7430	TheraPEAK™ X-VIVO™ 10 Medium With L-glutamine, Without Gentamicin or Phenol Red, With Human Transferrin	1 L		
BEBP02-055Q	BEBP02-055Q	TheraPEAK™ X-VIVO™ 10 Medium With L-glutamine, Without Gentamicin or Phenol Red, With Recombinant Transferrin	1 L		
	BE02-060F*	X-VIVO™ 15 Medium With L-glutamine, Gentamicin and Phenol Red, With Human Transferrin	500 mL		
04-4180	BE02-060Q	X-VIVO™ 15 Medium With L-glutamine, Gentamicin and Phenol Red, With Human Transferrin	1 L		
BE02-053Q	BE02-053Q	X-VIVO™ 15 Medium With L-glutamine, Gentamicin and Phenol Red, With Recombinant Transferrin			
BP04-744Q	BEBP02-061Q	heraPEAK™ X-VIVO™ 15 Medium With L-glutamine, Without Gentamicin or Phenol Red, With Human Transferrin			
08-879H	BE08-879H	TheraPEAK X-VIVO™ 15 Medium with L-Glutamine, without Gentamicin and Phenol Red, 5L bag	5 L		
BEBP02-054Q	BEBP02-054Q	TheraPEAK™ X-VIVO™ 15 Medium With L-glutamine, Without Gentamicin or Phenol Red, With Recombinant Transferrin			
04-4480	BE04-448Q	X-VIVO™ 20 Chemically Defined, Serum-free Hematopoietic Cell Medium, With L-glutamine, Gentamicin and Phenol Red	1 L		

^{*} Europe only

CHO Media

PowerCHO™ Chemically Defined, Serum-free CHO Media

- Chemically defined PowerCHO™ Media bring new levels of cell proliferation and protein production to chemically defined media
- Maintain high viability (>90%) at high cell densities


ProCHO™ Protein-free CHO Media

- Multiple formulas to optimize your protein-free applications including adherent and suspension cells, with high proliferation rates and high protein yield
- Directly convert cultures from adherent with serum to suspension without serum


PowerCHO™ Advance Medium

- Chemically defined, no raw materials of animal origin, serum-free and manufactured to regulatory standards
- Designed for growing and feeding CHO cells in serum-free conditions
- Allow for easier filtration while maintaining cell growth and viability.
- Provides protein titers equivalent or better compared to competitors

PowerCHO™ and ProCHO™ IgG Production



PowerCHO™ Advance Productivity Compared to Competitive Media



	General use	Non-animal origin	Protein-free	Chemically defined	Optimized for
PowerCHO™ 1, 2, CD Media	_				Suspension CHO
ProCHO™ 4, 5 Media	-			_	Suspension CHO
PowerCH0™ Advance Medium	_		•	•	Suspension CHO

Ordering information

Cat. No US	Cat.No EU	Product	Size
BP12-7700	BP12-770Q	PowerCHO™ 1 Chemically Defined Medium	1 L
BP12-771Q	BELN12-7710	PowerCHO™ 2 Chemically Defined Medium	1 L
12-771P10	BE12-771P10	PowerCHO™ 2 Chemically Defined Medium	10 L
12-771P20	BE12-771P20	PowerCHO™ 2 Chemically Defined Medium	20 L
BEBP12-029Q	BEBP12-029Q	ProCH0™ 4 Medium	1 L
BP04-919Q	BP04-9190	ProCHO™ 4 Medium, without Phenol Red	1 L
BP12-766Q	BELN12-7660	ProCHO™ 5 Medium	1 L
BE12-766P10	BE12-766P10	ProCHO™ 5 Medium	10 L
BE12-766P20	BE12-766P20	ProCHO™ 5 Medium	20 L
BE12-927Q	BE12-927Q	PowerCHO™ Advance Medium	1 L
02-0440	BE02-044Q	PowerFeed A	1 L
BE02-052Q	BE02-052Q	PowerFeed A with lipids	1 L
BE02-056Q	BE02-056Q	CH0 Xtreme™ Feed, Chemically Defined	1 L
BEBP17-855E	BEBP17-855E	ProHT™ Supplement 100x	100 mL

ProCHO™, PowerCHO™ Media and their Feeds are also available in powder format, contact Lonza for information.

Renal Media

Pro293™ Chemically Defined Media

- Formulas for both suspension and adherent culture

ProVero™ 1 Medium

- Supports Vero and MDCK in viral vaccine applications


UltraMDCK™ Chemically Defined Medium

 Supports MDCK in viral vaccine applications at both high and low densities


ProMDCK™ Medium

 Serum-free medium that supports the expansion and virus infection of MDCK cells in planar culture and on microcarriers

Pro293™s Easy Adaptation to Suspension Culture



Pro293™a Supports Adherent High-density Cell Growth



	General use	Non-animal origin	Protein-free	Chemically defined	Optimized for
Pro293™s Chemically Defined Medium	-				HEK 293 suspension
Pro293™a Chemically Defined Medium					HEK 293 adherent
ProVero™ 1 Medium				_	Vero, MDCK
ProMDCK™ Medium				_	MDCK in 2D and 3D
UltraMDCK™ Chemically Defined Medium					MDCK, Vero

Ordering information

Cat. No US	Cat. No EU	Product	Size
BP12-765Q	BEBP02-025Q	Pro293™s Chemically Defined Medium	1 L
BEBP12-764Q	BEBP12-764Q	Pro293™a Chemically Defined Medium	1L
BEBP02-030Q	BEBP02-0300	ProVero™ 1 Medium	1L
BE12-924Q	BE12-924Q	ProMDCK™ 2D Medium	1L
BE12-925Q	BE12-925Q	ProMDCK™ 3D Medium	1L
BEBP12-749Q	BEBP12-749Q	UltraMDCK™ Chemically Defined Medium	1 L

E

Insect Medium

Insect-XPRESS™ Protein-free Insect Cell Medium

 Versatile medium for shaker flask or stationary culture of SF9, SF21, High Five™ and Drosophila cells

	General use	Non-animal origin	Protein-free	Chemically defined	Optimized for
Insect-XPRESS™ Medium	_	-		-	SF9, SF21, High Five™ and Drosophila cells

Ordering information

Cat. No US	Cat. No EU	Product	Size
BP12-730Q	BELN12-730Q	Insect-XPRESS™ Medium	1 L
BEBP12-730P10	BEBP12-730P10	Insect-XPRESS™ Medium	10 L

Hybridoma Media

UltraDOMA-PF™ Chemically Defined, Protein-free Hybridoma Medium

 Based on the original hybridoma media for murine, human, and chimeric origin cells in a completely defined medium with no peptides or tissue extracts

ProDoma™ Serum-free Hybridoma Medium

- Non-animal, protein-free medium for your hybridoma needs
- Scalable from small flasks to large-scale bioreactors

	General use	Non-animal origin	Protein-free	Chemically defined	Optimized for
UltraD0MA-PF™ Chemically Defined Medium	•			-	Murine, human, and chimeric cell lines
ProDoma™ 3 Chemically Defined Medium	_			_	Murine, rat, human, and chimeric cell lines

Ordering information

Cat. No US	Cat. No EU	Product	Size
BP12-727F	BP12-727F	UltraDOMA-PF™ Chemically Defined Medium	500 mL
BEBP02-032Q	BEBP02-0320	ProDoma™ 3 Chemically Defined Medium	1 L

NSO Medium

ProNSO™ Chemically Defined, Protein-free Medium

- Designed for cultivation of NSO myeloma cells
- Maximise protein production by titration of ProNSO™ Lipid
 Chemically Defined Supplement

Ordering information

Cat. No US	Cat. No EU	Product	Size
BP12-7730	BP12-7730	ProNSO™ 1 Protein-free Medium, Chemically Defined	1 L
BP12-775J	BP12-775J	ProNSO™ Lipid Supplement - Chemically Defined	5 mL

Contact Information

North America

Customer Service: +1 800 638 8174 (toll free) Fax: +1 301 845 8338

order.us@lonza.com

Scientific Support: +1 800 521 0390 (toll free)

scientific.support@lonza.com

Europe

Customer Service: +32 87 321 611

order.europe@lonza.com

Scientific Support: +32 87 321 611 scientific.support.eu@lonza.com

International

Contact your local Lonza distributor Customer Service: +1 301 898 7025 scientific.support@lonza.com

International Offices

+61 3 9550 0883 Australia Belgium +32 87 321 611 Brazil +55 11 2069 8800 0800 91 19 81 (toll free) France Germany 0800 182 52 87 (toll free) India +91 22 4342 4000 Japan +81 3 6264 0660 +32 87 321 611 Luxemburg +65 6521 4379 Singapore

The Netherlands 0800 022 4525 (toll free)
United Kingdom 0808 234 97 88 (toll free)

Lonza Cologne GmbH - 50829 Cologne Germany

For research use only. Not for use in diagnostic procedures. PER.C6 and Permexcis are trademarks of Crucell N.V.

High Five is a trademark of Boyce Thompson Institute for Plant Research Inc.

All trademarks belong to Lonza or its affiliates or to their respective third party owners. The information contained herein is believed to be correct and corresponds to the latest state of scientific and technical knowledge. However, no warranty is made, either expressed or implied, regarding its accuracy or the results to be obtained from the use of such information and no warranty is expressed or implied concerning the use of these products. The buyer assumes all risks of use and/or handling. Any user must make his own determination and satisfy himself that the products supplied by Lonza Group Ltd or its affiliates and the information and recommendations given by Lonza Group Ltd or its affiliates are (i) suitable for intended process or purpose, (ii) in compliance with environmental, health and safety regulations, and (iii) will not infringe any third party's intellectual property rights.

©2019 Lonza. All rights reserved. CD-BR020 07/19

www.lonza.com/specialty-media