

Datasheet

Catalog Number: NBAB-160170

Version: RN6.0



Product Name	NebuSelect™ Anti-Protein G Polyclonal Antibody
Alias/Synonyms	Protein G
Description	NebuSelect™ Anti-Protein G Polyclonal Antibody (CAT# NBAB-160170) was produced in rabbits immunized with Native Full Length Protein corresponding to Streptococcus sp. 'group G' spg, and it reacts with Bacteria. This antibody is validated for ELISA, IF/ICC, WB.
Host / Isotype	Rabbit / IgG
Type	Polyclonal
Immunogen	Recombinant Protein G
Conjugate	Unconjugated
Tested Applications	ELISA, IF/ICC, WB
Dilution	ELISA 1:5000 ~ 1:20000 IF/ICC 1:50 ~ 1:500 WB 1:5000 ~ 1:20000
Tested Reactivity	Protein G
Buffer	0.2 µm filtered solution in TBS, pH7.4, containing 0.02% NaN3
Size	100ul; 500ul; 1ml
Storage	<p>This liquid antibody can be stored at 2°C-8°C for at least two weeks. Avoid repeated freeze/thaw cycles.</p> <p>The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or antibody diluent, the antibody is stable for at least two weeks at 2-4°C.</p>